

Comfort Shield



Comfort Shields were developed and deployed to protect our nurses and doctors during the COVID-19 pandemic. As we continue to support and provide protection to frontline workers all over the world, we want to expand protection to businesses preparing to reopen. Restrictions will be lifted, but businesses will want to provide PPE equipment to their employees.

Dentists, chiropractors, therapists, restaurant & bar workers, ride share drivers, gas station attendants, bank tellers, grocery store clerks, teachers, retail store associates... The list of those in need of the Comfort Shield is endless. Protect your people with an American made product.

MATERIAL

- PETG .020" Clear Plastic Face Shield
- FDA medical grade headband

CLEANING INSTRUCTIONS

- Use alcohol to clean the blue headband
- Use disinfectant spray with a soft cloth on the plastic face shield

SHIPPING SPECS

- Comfort Shields ship flat
- 25 Comfort Shields fit in a 14"x11"x5.5" box, weighing only 5.5 lb.

30-SECOND ASSEMBLY

- Match up the top 4 holes on the clear shield with the 4 prongs on the headband
- Snap the 4 prongs into the 4 holes on the shield
- After you have successfully attached the headband to the shield, remove the protective adhesive film



Angle View



Top View



Front View